# Home Oxygen Equipment & User Guide

Oxygen safety and  
maintenance of your  
equipment

## Welcome to Air Liquide Healthcare Home Oxygen Therapy


Air Liquide Healthcare is a world leader in medical gases, home healthcare and healthcare specialty products. Our aim is to provide our patients the continuum of care from hospital to home, we put our patients at the 'centre' of their healthcare, while treating you with dignity, respect and compassion.

*We are glad that we can partner with you in providing your Home Oxygen therapy needs.*

Our Home Oxygen program is designed to provide solutions to meet your individual requirements, as well as ongoing support to ensure that your Medical Oxygen therapy needs are being met today and over time. We are very happy to be able to help you be mobile and get the most out of life.

To download user manuals, instructions and brochures for our products all in one place visit the Product Manuals page on our website at any time that is convenient for you. At Air Liquide Healthcare we strongly believe that a high level of support is the key in ensuring the effectiveness of your therapy. We supply a range of Medical Oxygen supplies and accessories from Cannula, Comfort Solutions, Connectors & Extension Tubing, which can be ordered conveniently online via or eShop.

Our services include

Thank you again for choosing to partner with us for your medical Oxygen needs. We look forward to assisting you in getting the most out of each breath.

Rely on world-class expertise to assist you in getting the most out of each breath.

Our experienced nationwide team of medical-technical experts and knowledgeable customer service staff ensure peace of mind for all home Oxygen users, aged care facilities, area health and government agencies.

**Trust over 25 years of proven reliability.**

*Follow us* @airliquidehealthcareau on the below channels!

## How do I contact Air Liquide Healthcare?

✉ [alhenquiries@airliquide.com](mailto:alhenquiries@airliquide.com)

☎ 1300 360 202

🌐 [au.healthcare.airliquide.com/oxygen](http://au.healthcare.airliquide.com/oxygen)

🛒 [store.airliquidehealthcare.com.au](http://store.airliquidehealthcare.com.au)

# CONTENTS

<b>General Oxygen Safety</b>	4
<b>What equipment is used to supply my home Oxygen therapy?</b>	6
<b>How do I travel with my Oxygen equipment?</b>	10
<b>Home Oxygen Products</b>	11
<b>Regulators and Oxygen Conserving Devices</b>	13
<b>Respiratory Support Products</b>	14
<b>Stationary Oxygen Concentrators</b>	15
- Stationary Oxygen Concentrator safety	
- Using your stationary Oxygen concentrator	
- Cleaning your stationary concentrator	
- Concentrator alarms	
- Fixing problems	
<b>Portable Oxygen Concentrators</b>	20
- POC safety	
- Fixing problems	
<b>Medical Oxygen Bottles</b>	23
- Oxygen Bottle safety	
- Using your Oxygen bottle	
- How long will your Oxygen bottle last	
- Going out with your Oxygen bottles	
- Fixing problems	
<b>Paediatric Equipment</b>	29
<b>Oxygen Accessories</b>	31
<b>General Oxygen Equipment Cleaning</b>	33
<b>Therapy Installation and Subsequent Visits to a Patient's Home</b>	34
<b>Oxygen Reference Guide</b>	39
<b>Power failure emergency plan</b>	45
<b>My emergency plan</b>	46

# General Oxygen safety

Before operating your Oxygen equipment, make sure you read the safety information in this booklet. It will help you to operate your equipment safely and answer many questions.

Only use your medical Oxygen at the flow or setting prescribed by your doctor when at rest, exercising and while sleeping.

- Fire needs Oxygen to burn

Oxygen assists in making things burn much faster


These safety rules apply whenever Oxygen is in use  
Please keep this booklet handy for your reference.

Read this information prior to using or helping others use their Oxygen equipment

For your safety please read the following before commencing your oxygen therapy

**NO SMOKING** when oxygen is in use This includes electronic cigarettes

**DO NOT** use oxygen within three metres of any open flame

- NO lit candles
- NO gas pilot lights
- NO gas or wood burning stoves
- NO open fireplaces


Put 'OXYGEN IN USE' sign where you use and keep your oxygen. Make sure others can see the sign

Put the sign on the:

- Front door of your house
- Garage door
- Car window

Air Liquide Healthcare will give you as many signs as you need

Clothes and blankets can easily 'absorb' oxygen, making them burn easier and faster  
Use oxygen in rooms where there is lots of fresh air

- **DO NOT** cover your oxygen equipment with blankets, towels or clothing
- **DO NOT** store your oxygen in a small or cluttered area
- **DO NOT** store outside or in an unsecured area

Turn your oxygen **OFF** when you are not using it


Make sure your home has at least one smoke alarm Check it's batteries are replaced regularly!

NO SPARKS

Check that ALL of your electrical appliances are in good working order

**DO NOT** use electrical equipment that get hot or spark near your oxygen

Use your oxygen equipment in the upright position

Make sure your oxygen system does not fall over by ensuring it is adequately secured

For safety, cook with a microwave

If you use a stove, turn your oxygen OFF while you cook

Make sure your oxygen system does not get wet

- **DO NOT** store your oxygen in moist or humid areas
- **DO NOT** store your oxygen in the bathroom
- **DO NOT** store your oxygen in the laundry room
- **DO NOT** store your oxygen outside


If your nose is dry, Air Liquide Healthcare or your doctor can suggest special creams

- **DO NOT** use grease, oil, oil-based or petroleum-based skin creams (e.g. vaseline, paw paw ointments, gels, etc) when using oxygen
- Look at the "Oxygen Accessories" section of this booklet for compatible products you can use

Only use oxygen tubing provided to you by Air Liquide Healthcare  
Long oxygen tubing can help you get around your home  
Be careful not to trip

Only clean with a cloth

Keep your oxygen system clean and dust free using a damp cloth

**DO NOT WASH** your oxygen system directly with water

# WHAT EQUIPMENT IS USED TO SUPPLY HOME OXYGEN THERAPY?

Your physician will decide whether you require Oxygen at home and/or when you leave the house. Air Liquide Healthcare will help you select the equipment that best suits your needs. Oxygen can be provided through three main supply systems. These systems may be used alone or in different combinations.

Systems include:

1. stationary Oxygen concentrators
2. portable Oxygen concentrators
3. Oxygen cylinders

## WHAT IS A STATIONARY OXYGEN CONCENTRATOR?

An Oxygen concentrator is an electrically driven device that extracts the Oxygen from room air. This means that the concentrator never runs out of Oxygen and never needs to be refilled. Most stationary Oxygen concentrators weigh anywhere between 18 – 30 kg and cannot be used outside of the home.

Stationary Oxygen concentrators feature Oxygen sensing devices that will warn you if the purity of the Oxygen delivered decreases, so that you can use an alternative Oxygen supply system and contact Air Liquide Healthcare for assistance.

## WHAT IS A PORTABLE OXYGEN CONCENTRATOR?

Portable Oxygen concentrators (POC) operate in a similar way to stationary Oxygen concentrators but are designed to be used while you are outside of the home. The weight of portable Oxygen concentrators varies between models, ranging from 2 - 8kg. As a general rule, heavier pocs are capable of delivering higher Oxygen flow rates.

The battery duration for portable Oxygen concentrators varies between models. Batteries are recharged from your mains power at home, and most have connections to recharge from your vehicle. Your poc should be serviced regularly to ensure that it is working properly and that the Oxygen purity stays within the manufacturer's recommendations.

**Since the amount of Oxygen produced varies greatly between models, you will need to be tested on the specific device you intend to use to ensure that it will meet your Oxygen needs.**

# WHAT ARE OXYGEN CYLINDERS?

Oxygen cylinders contain Oxygen in gas form. they are filled under pressure to contain a greater volume of gas. The pressure is released with the use of a regulator which is also used to set the prescribed flow rate.

Cylinders can also be used with an Oxygen conserving device which only delivers a flow of Oxygen when you breathe in. Your cylinder will last much longer since there is no Oxygen flowing when you breathe out.

Oxygen cylinders are available in a variety of sizes and are separated into two main groups; portable and stationary.

Portable (small) cylinders are lightweight cylinders used outside of the home, and are available in a variety of sizes. Smaller models are contained in a bag which can be carried over your shoulder or used as a backpack.

The larger models can be transported using a trolley. Stationary (large) cylinders are normally used as a safety 'back up' for stationary Oxygen concentrators in the event of power failure or equipment malfunction. They are available in a variety of sizes. Selection is based on your prescribed flow rate and how far away you live from an Air Liquide Healthcare office.

**Your physician will provide your medical Oxygen prescription.**

**Air Liquide Healthcare will assist you in selecting the equipment best suited to your needs.**

# WHAT ACCESSORIES DO I NEED WITH MY OXYGEN THERAPY?

Oxygen is typically delivered from the supply system to your nose with a 'cannula'. Extension tubes are available in a variety of lengths (up to 11 metres) to allow you to move about the house while connected to your Oxygen supply system. Both the cannula and the extension tubing are available in a variety of sizes and lengths. Based on your prescription, Air Liquide Healthcare will assist you in choosing accessories.

Firesafe cannula valves\* can be connected to your Oxygen accessories to provide an extra line of defence, that can limit the impact of an Oxygen fire by automatically arresting the Oxygen flow in the event of a fire.

Other specialty accessories are available as prescribed by your physician.

# HOW OFTEN DO I REPLACE MY OXYGEN ACCESSORIES?

Your cannula will become dirty and brittle with use and may hurt your ears and skin if not replaced at least once to twice a month.

Extension tubing is disposable and should be replaced every three to six months to ensure that no leaks occur.

The firesafe valve is typically replaced every six months.

# HOW DO I USED OXYGEN IN MY EVERYDAY LIFE?

Integrating Oxygen therapy into your daily routine will not only prolong your life but enhance your quality of life(1)(2). Remember that the objective of using Oxygen therapy is to increase the amount of Oxygen in your blood, supplying enough to all organs during rest and during physical activities.

Use your Oxygen as prescribed while you go about your normal routine such as running errands, doing housework, watching television or reading a book. It is important to use your Oxygen even when outside of the home. Use Oxygen as prescribed in your everyday life, when:

## WHAT CAN I DO IF MY EARS OR NOSE ARE SORE?

If your nose becomes dry or sore you can use a product such as Nozoil, which is safe to use with Oxygen therapy.

You should never use any petroleum based products (such as creams and oils) with Oxygen.

If your ears are sore, you may want to try a product such as EZ Wraps soft foam tubes that fit over the cannula. They provide a soft cushioning for the sensitive tissue around your ears.

If you are experiencing discomfort from your Oxygen therapy contact Air Liquide Healthcare for assistance.

## ARE THERE ANY SUPPORT GROUPS FOR PEOPLE ON OXYGEN THERAPY?

Support groups are invaluable for Oxygen users, their carers and family. Meeting people going through similar life changes can be encouraging and enlightening. Support groups offer an opportunity to network, learn and share skills to help in a range of issues associated with health problems such as feeling low, irritability, isolation and confusion. There is also a strong focus on how to stay healthy and active.

There are a variety of established support groups including support from your area health service and Lung foundation Australia. For information on a lung support group nearest you, please contact your local healthcare facility or Lung foundation Australia:

<http://lungfoundation.com.au/patient-area/patient-support>

# HOW DO I TRAVEL WITH MY OXYGEN EQUIPMENT?

Air Liquide Healthcare can provide Oxygen to most popular destinations and will advise you of any associated costs. Always try and plan the arrival to your destination on a week day, and follow the pre-travel planning check list below.

## YOUR PRE-TRAVEL PLANNING CHECKLIST

Planning well in advance is the most important part of ensuring a successful trip. Your travel plan should include all parties involved in making arrangements for the supply of your Oxygen so that you can travel safely to and from your destination.

Before you travel ensure:

- You have a current copy of your Oxygen prescription.
- You have consulted with your airline, train and/or cruise liner regarding any special medical clearance forms required to be completed by your physician before departure. Allow at least 4 – 6 weeks for this process.
- You are aware of the total time of your trip so that you can ensure adequate supply of Oxygen.
- You have discussed with your physician the use of mobility devices such as portable Oxygen concentrators or cylinders with Oxygen conserving devices.
- You have the contact details of the Air Liquide Healthcare outlet.
- You have spare accessories such as nasal cannulas/masks, tubing and connectors.

## KEY POINTS

### **Notify Air Liquide Healthcare at least:**

- two weeks in advance for destinations within Australia.
- one month in advance for overseas destinations.

### **Air Liquide Healthcare will require the following information:**

- travel dates and mode of travel (train, plane etc.).
- Length/time of journey.
- travel details, including any flight, cruise or train references.
- travel itinerary, including all addresses of hotels and travel destinations.
- details of a contact person during travel.

# HOME OXYGEN PRODUCTS

Air Liquide Healthcare has a wide range of Oxygen equipment available for rental and purchase. For more product information and product manuals visit our website at <http://home-oxygen.com.au/products/>

## STATIONARY OXYGEN CONCENTRATORS

These devices are electrically powered and supply Oxygen extracted from room air so they never run out. All of Air Liquide Healthcare's Oxygen concentrators feature built-in Oxygen monitors which alarm in the event of low Oxygen purity and power disconnection.

## PORTABLE OXYGEN CONCENTRATORS (POCS)

Unlike stationary Oxygen concentrators which deliver a continuous flow of Oxygen, pocs use pulse dose technology which only provides Oxygen when you breathe in. Some models may also feature a continuous flow option. Pocs are lightweight and are ideal to help you get the benefits of Oxygen while you are outside of the home.

Air Liquide Healthcare has a range of portable and transportable Oxygen concentrators available:

- Portable Oxygen concentrators (poc) are devices which can be easily carried or moved and typically weigh less than 2.5 kg.
- Transportable Oxygen concentrators (toc) are devices that normally weigh over 2.5 kg and utilise a trolley for transportation. the advantage of these devices is that they tend to have a greater Oxygen output when compared to portable devices.

**Please note that portable/transportable Oxygen concentrators are not intended to replace a stationary Oxygen concentrator. Most manufacturers recommend that a stationary Oxygen concentrator is used at home, with pocs/tocs only being used outside of the home.**

# OXYGEN CYLINDERS

## PORTABLE OXYGEN CYLINDERS

Air Liquide Healthcare has an extensive range of portable Oxygen cylinders available for rent. Cylinders vary in size and in the volume of Oxygen they contain. our popular 'b' size (160L) cylinder is the smallest and lightest Oxygen cylinder available in Australia. This cylinder is very popular with people who will benefit from using an Oxygen cylinder while they are mobile.


Portable Oxygen cylinders are contained in a robust carry bag or can be transported using a trolley depending on the cylinder size/weight.

## STATIONARY OXYGEN CYLINDERS

Stationary Oxygen cylinders are recommended as an emergency backup source of Oxygen in case of power failure or Oxygen concentrator malfunction. As a general rule 'D', PRC and TAKE o2® 5.0 size cylinders are provided in urban areas. 'E', PRE and TAKE o2® 20 size cylinders are provided in regional areas. Stationary Oxygen cylinders are also available with integrated regulators and flow meters.

## OXYGEN CYLINDER GUIDE

# REGULATORS AND OXYGEN CONSERVING DEVICES

## OXYGEN REGULATORS

Regulators are pressure-regulating devices that attach to Oxygen cylinders. Prescribed flow rates are set with an adjustable flow meter. Regulators deliver a continuous flow of Oxygen and are available in low flow and high flow ranges

## OXYGEN CONSERVING DEVICES (OCD)

An Oxygen conserving device (ocd) is a battery operated regulator that is used to extend the duration of use of portable Oxygen cylinders.

## OXYGEN CANNULAS AND MASKS

Air Liquide Healthcare carries a variety of Oxygen cannula and tubing to suit all of your Oxygen delivery needs. Our goal is to ensure the most comfortable and effective means of Oxygen delivery throughout your daily activities.

The use of an ocd with a portable Oxygen cylinder can increase the duration of a cylinder by up to 7 times, providing economic benefits as well as increased mobility. OCDs use pulse-dose technology. This technology delivers a precise pre-determined volume of Oxygen during the inhalation stage only of the breathing cycle. Some OCDs can be adjusted to continuous flow mode if necessary (for e.g. if the battery runs out).

## COMFORT PRODUCTS

We also provide a range of products to improve the comfort of your Oxygen therapy, including: Ear cushions, which help relieve soreness around your ears caused by a Oxygen cannula and mask. As well as Nozoil™ - this all natural oil relieves symptoms of nasal dryness and is compatible for use with Oxygen therapy.

## FIRESAFE™ VALVES

Firesafe™ valves can be connected between your Oxygen accessories and the Oxygen source to provide an extra line of defense. Firesafe™ valves can limit the impact of an Oxygen fire by automatically arresting the Oxygen flow in the event of a fire.

# RESPIRATORY SUPPORT PRODUCTS

Air Liquide Healthcare carries a range of products to assist patients with COPD and Asthma, including:

## AEROBIKA™

The AEROBIKA™ oscillating positive expiratory pressure (opep) therapy system works with the lungs to help loosen and remove mucus. The AEROBIKA™ may be used as recommended by health professionals for people diagnosed with respiratory or other recurring conditions which may affect the lungs.

The AEROBIKA™ is the only opep device that can be used simultaneously with the AEROECLIPSE™ breath actuated medication nebuliser, resulting in reduced treatment times. The AEROBIKA™ is non positional.

Always check with your healthcare professional to make sure the AEROBIKA™ is suitable for you.

## Stationary aerosol compressors

Suitable to nebulise all respiratory medications.

## Aerosol spacers

Used to delivery respiratory medication with most 'puffers', or metered dose inhalers (mdi's)

## AEROECLIPSE™

The AEROECLIPSE™ BAN nebulizer is a breath-actuated small volume medication nebulizer. Aerosolized medication is only produced when you inhale through the device. This means between breaths or during breaks in treatment, prescribed medication is contained in the cup until you are ready to inhale it. The AEROECLIPSE™ BAN has a unique patient friendly green feedback button which moves down on inhalation. The AEROECLIPSE™ BAN can be used simultaneously with the AEROBIKA™. Reducing therapy time may improve treatment compliance.

## Peak flow meters

Recommended for adults and children over the age of 6:

- for people with severe airway limitations.
- for recording the effect on your lungs when your doctor changes your medication.

## Nebuliser medication kits

May be used with most aerosol compressors

# SUCTION EQUIPMENT

## SUCTION COMPRESSOR

The Aspira Go Suction compressors are available portable models.

The portable units operate from mains power or lithium battery, making them ideal for those who require assistance outside of the home and in an emergency. The portable compressors are equipped with a vacuum regulation knob and a vacuum gauge. The collection jar has a protection device preventing entrance of liquid in the suction pump through a float shut-off located in the collection jar lid

## SUCTION ACCESSORIES

### Suction Catheters

Suction catheters feature transparent medical grade plastic with smooth grooved suction tip and regulator control vents on the handle making them safe and easy to use. A large selection of catheter sizes are available.

### Yankauer Suction

This catheter combines a sterile, single use tip with sturdy shatter resistant construction to enable retraction use. Large lumen resists clogging and the transparency allows clear visualisation of liquids.

**\*Please note that not all equipment listed in this document is included in your funding. Please check with your funding body or call Air Liquide Healthcare**

# STATIONARY OXYGEN CONCENTRATORS

## Stationary Oxygen Concentrator Safety

Your home healthcare technician will show you how to use your oxygen concentrator safely

- Make sure you ask lots of questions
- You will be asked to sign a check sheet to say you understand how to use your Oxygen therapy equipment safely
- Others helping you with your Oxygen concentrator should be there when your Oxygen is being set up

Your Oxygen concentrator should be plugged directly into a wall outlet

- **DO NOT** use electrical extension cords
- **DO NOT** use power boards or double adapters

- A FireSafe™ Valve is a safety device that stops the flow of Oxygen when exposed to fire. It is used with Oxygen concentrators
- Some Oxygen concentrator models have built in FireSafe™ Valves
- If your concentrator does not have a FireSafe™ Valve, Air Liquide Healthcare will provide you one during setup. The Valve will need to be connected to a short tube where the Oxygen comes out. Your Oxygen tubing is then connected to the FireSafe™ Valve

**DO NOT COVER**

Keep all sides of your Oxygen concentrator at least 30 cm away from:

- Furniture
- Walls
- Curtains

Take care not to block the filter.

**DO NOT** cover your concentrator with blankets, towels or clothes

Always keep your concentrator on a flat surface **DO NOT** place your Oxygen concentrator on or near a staircase

Oxygen concentrators are electrical appliances, so if you cannot be without Oxygen, you should have a backup Oxygen system

**DO NOT** use your concentrator when it is colder than 5°C or warmer than 35°C in the house

# Using your stationary Oxygen concentrator

**1** Plug your Oxygen concentrator directly into the wall outlet


**2** Turn the switch on the front to 'ON'  
An alarm will sound every time you turn your concentrator on. This can last up to 60 seconds  
Once powered the alarm will stop. If there is no electrical power and the concentrator is turned on, the concentrator will not stop alarming  
The alarm is quite loud so that you can hear it from around the house (or if you are sleeping)

**3** Set your Oxygen flow  
Your doctor will have prescribed this in litres per minute (L/min)

- Turn the knob 'left' to increase the flow
- Turn the knob 'right' to decrease the flow

Read the ball at eye level so that the ball is in the centre of the line  
DO NOT change the flow without instructions from your doctor

**4** Connect your Oxygen tubing  
Connect your long Oxygen tubing to the outlet of the concentrator. DO NOT use more than 15m of long tubing  
If using a FireSafe™ Valve, connect the long tubing to the valve  
Connect the swivel to your long tubing  
Connect the cannula to the swivel  
Your cannula should be worn with the curved ends pointing down. The tubing goes around your ears and under your chin. Tighten it slightly but make sure you are comfortable

**5** Turn the Oxygen 'OFF' when you are not using it  
When you turn your Oxygen concentrator off and turn it back on, the flow setting will go back to the last setting

## Cleaning your stationary Oxygen concentrator

- If your Oxygen concentrator has a filter, it must be cleaned every week
- The filter is normally found on the side or back of the concentrator
- Rinse the filter under warm soapy water and pat it dry with a towel
- The filter must be COMPLETELY DRY before you put it back in the machine
- **DO NOT** turn your machine on without the filter in place

Keep your concentrator dust free by wiping it weekly with a clean damp cloth

Never open the cabinet of a concentrator

## Moving your stationary Oxygen concentrator

**Oxygen concentrators are heavy. Always get help to move them**

Always transport your Oxygen concentrator in an upright position

If you are travelling by car, have a friend or family member lift it onto the back seat

Secure it upright to the backseat of your car with a seatbelt

**DO NOT** put your stationary Oxygen concentrator in the boot of your vehicle

## Concentrator alarms

The Oxygen concentrator alarm will flash and make a loud beeping sound.

### **The alarm will sound if:**

- The concentrator is turned on and no electrical power is available
- The Oxygen concentrator is not making enough Oxygen If you hear it alarming you should:

- Check that your Oxygen concentrator is plugged into the wall
- Check for power around your home

### **If it is not a power issue, you should:**

- Turn the Oxygen concentrator off
- Use your back up Oxygen option (if you have one)
- Call Air Liquide Healthcare for assistance

# Fixing problems

## Problem 1: Alarm sounds

Possible cause	Solution
Low Oxygen concentration - light on	<ul style="list-style-type: none"><li>• Check that the Oxygen concentrator is getting electrical power. If you have power to the concentrator, the alarm may mean the Oxygen concentrator is not making enough Oxygen</li><li>• Use your back up Oxygen option if you have one and call Air Liquide Healthcare</li></ul>
Power off	Check that the connection at your wall switch is on
Circuit breaker tripped or fuse blown	Check your house fuse or reset your circuit breaker
Internal problem with the unit	Use your back up Oxygen source if you have one and call Air Liquide Healthcare

## Problem 2: No Oxygen flow from the concentrator - No humidifier in use

Possible cause	Solution
Tubing is disconnected	<ul style="list-style-type: none"><li>• Connect the tubing</li><li>• Check that the ball is going to the correct setting on the flow meter and adjust if needed</li><li>• Check the tubing starting from the concentrator to your nose</li></ul>
The tubing is damaged or leaking	Replace it with new tubing
The flow meter is turned off	Turn the flow meter on and check that your machine is set to your prescribed flow
The Oxygen concentrator is turned off	Turn the machine on and check the main power at your wall switch

### Problem 3: No Oxygen flow from the concentrator - With humidifier in use

Problem cause	Solution
The humidifier is leaking	Tighten the humidifier lid. Take care not to over tighten
The humidifier jar is cross-threaded on the lid	Disconnect the humidifier jar from the lid and reattach. Take care not to over tighten it
The humidifier jar is cracked	Replace the humidifier
You hear a 'whistling' noise from your tubing	<p>Drain any water in the tubing</p> <p>Consider using a water trap to collect condensation</p> <p>Look at the "Oxygen Accessories" section of this booklet for more information on water traps</p>
The humidifier is making a 'popping' noise	Your tubing may be blocked or kinked. Check that no heavy objects are on your tubing

**DO NOT** open the Oxygen concentrator  
**DO NOT** try to fix your Oxygen concentrator on your own  
 Use your back up Oxygen source if you have one.  
 Call Air Liquide Healthcare for assistance

# PORTABLE OXYGEN CONCENTRATORS (POC)

## POC safety

The same general Oxygen safety and stationary Oxygen concentrator safety rules apply to POCs. Below are additional precautions you should take to stay safe

Be careful not to get your POC wet  
DO NOT use your POC in the rain or when bathing

Most POCs cannot be used with extended tubing.  
Check your POC user guide

only use nasal cannulas with your POC  
Oxygen masks and special accessories cannot be used with POCs

DO NOT use a humidifier with your POC

Charge your POC batteries every month even if you don't use them

If your POC has an air inlet filter, always have a spare one

- **DO NOT** use your POC without the filter in place
- **DO NOT** block your filter or vent

When using your POC in the car, always turn the car on before plugging it into the cigarette lighter

If you have a new POC, power up the battery of your POC completely before using it for the first time  
This can take up to half a day

# Fixing problems

## Problem 1: It won't turn on

Possible cause	Solution
The battery is flat or poorly charged	Charge your battery
The battery has not been fitted correctly	Take the battery out Check your user guide and put it back in
POC may have gotten too hot or too cold	Let your POC reach a normal temperature and try again
Bad connection with your car cigarette lighter	Check that the DC cord is a proper fit for your car. Try plugging it back in

## Problem 2: The POC won't deliver a 'puff'

Possible cause	Solution
The POC is not turned on	Turn the POC on
Tubing is bent, twisted or too long	Check your tubing to make sure it is not bent or twisted. 1. Straighten out the tubing 2. Check if it is connected properly to the POC 3. Replace the tubing if required 4. Check your user guide to make sure your tubing is not too long
The flow meter is turned off	You may be breathing too quickly. Slow your breathing down and make sure you are breathing through your nose

### Problem 3: Alarm sounds

Possible cause	Solution
The POC is not ready to make Oxygen	If you have just turned the POC on, wait a few minutes for it to warm up and start making Oxygen
Technical fault with the POC	Restart your POC Try and replace or recharge your battery If this does not work call Air Liquide Healthcare
Low battery	Replace and/or recharge your battery Connect to mains power supply
You are breathing too quickly	Try to relax and slow your breathing down. The machine resets itself when you are breathing slower
The POC cannot detect your breathing	1. Check the connection from the tubing to your POC 2. Make sure the tubing is not kinked, bent or too long Replace if necessary

**DO NOT** open the POC  
**DO NOT** try to fix your POC on your own  
If these solutions don't work, call Air Liquide Healthcare for assistance

For further information on the safe use, operation and maintenance of your portable Oxygen concentrator, please refer to the supplier's user manual supplied with the device.

Supplier user manuals for all devices can also be found on our website [home-oxygen.com.au](http://home-oxygen.com.au)

# MEDICAL OXYGEN BOTTLES

## Oxygen bottle safety

The same general Oxygen safety rules apply to Oxygen bottles. Below are additional precautions to help you stay safe

### Storing your oxygen bottles

Store your bottles

- In a safe and secure area where they cannot be tampered with or stolen
- Away from areas with oil or grease such as kitchens and garages

Keep your spare Oxygen bottles in a well ventilated area

- **DO NOT** store them in a cluttered space
- **DO NOT** store them near curtains
- **DO NOT** store them near naked flames such as kitchens or BBQ areas
- **DO NOT** store them near flammable or combustible items

### Handling

Wash your hands before handling your Oxygen bottles **DO NOT** use hand creams

Be careful with your Oxygen bottle:

- **DO NOT** pick them up by the valve
- **DO NOT** drop them
- **DO NOT** let them fall over
- **DO NOT** remove the valve protector (where fitted)


Always close the Oxygen bottle valves when they are not in use

Only use Air Liquide Healthcare carry bags and/or trolleys

# Using your bottle safety


Before you begin, check your Oxygen bottle. DO NOT use your Oxygen bottle if:

- The **seal** is broken or missing
- There is no **heat tag**
- The **product label** is damaged, missing or incorrect
- It has **dents** or **scratches**

If you think something may be wrong with the Oxygen bottle do not use it. Keep it aside.

**DO NOT** turn your Oxygen on before the seal and valve plug is removed

- 1** Prepare your cylinder
  - Remove the seal and valve plug
  - Point your Oxygen bottle away from yourself or others
  - Turn the 'cap' at the top of your Oxygen bottle to the left to OPEN then CLOSE quickly
- 2** Putting a regulator or OCD onto your Oxygen bottle
  - Put the regulator or OCD (A) over your Oxygen bottle
  - Line up the pins on the regulator/OCD with the holes on the Oxygen bottle
  - Tighten the Regulator/OCD by turning knob (B) to the right. Do not over tighten
  - Slowly open the black 'cap' until it is completely OPEN
  - Set the Oxygen to the flow prescribed by your doctor using knob (C)

Always turn your Oxygen bottle off before you put your regulator or OCD on or take it off
- 3** Taking the Regulator/OCD off the Oxygen bottle
  - Turn to knob (C) to that the flow setting is at '0'
  - Turn the black 'cap' on the top of your Oxygen bottle to the right to CLOSE. The needle will slowly drop to '0' on the dial (D)
  - Loosen knob (B) by turning it to the left and lift the regulator/OCD off the Oxygen bottle to take it off

# How long will your Oxygen bottle last

How long your Oxygen bottle lasts depends on:

- The size of your Oxygen bottle
- If you are using a regulator or an OCD
- Your flow rate as prescribed by your doctor
- How fast you are breathing (Only if you are using an OCD)

Use the tables (on page 25) below to see how long your Oxygen bottle will last in HOURS.

Please note these are an estimate only.

Always check the dial on your regulator/OCD to make sure you have enough Oxygen left in your Oxygen bottle

### Table 1: Using an Oxygen Conserving Device (OCD) with Conventional Cylinders

OCDs give you Oxygen only when you breathe in. The times shown on the table are based on a breathing rate of 20 breaths per minute with the OCD set to 'pulse' mode.

'Please note that when the OCD is set to 'Continuous Flow (CF)' mode, it will deliver Oxygen at a fixed flow of 2 L/min. When in this mode use Table 2 to estimate bottle duration.'

OXYGEN BOTTLE SIZES <small>This is written on the product label of your Oxygen bottle</small>						
Flow Setting Pulse Dose	'B' 160L	'CS' 270L	'CH' 470L	'CL' 760L	'D' 1500L	'E' 4200L
1	10.2	17.3	30.1	48.7	96.1	276
2	5.8	9.8	17	27.5	54.3	156
3	4	6.8	11.9	19.2	37.9	108.5
4	3	5	8.7	14.1	27.8	79.6
5	2.4	4.1	7.2	11.5	22.7	65.1
6	2	3.3	5.7	9.3	18.4	52.7

### Table 2: Using a Paediatric Regulator (0 - 2LPM) with Conventional Cylinders

CYLINDER DURATION ESTIMATE IN HOURS AT VARYING FLOW SETTINGS						
Flow Setting with Continuous Flow underneath	'B' (160L)	'CS' (270L)	'CH' (470L)	'CL' (760L)	'D' (1500L)	'E' (4200L)
0.1l/pm	26	45	78	126	250	700
0.125l/pm	21	36	62	101	200	560
0.25l/pm	11	18	31	51	100	280
0.4l/pm	7	11	19	32	62	175
0.5l/pm	5	9	15	25	50	140
0.6l/pm	4	7.5	13	21	42	116
0.7l/pm	3.5	6.5	11	18	36	100
0.8l/pm	3	5.5	10	16	31	87
1.0l/pm	2.5	4.5	8	13	25	70
1.5l/pm	1.5	3	5	8	17	47
2.0l/pm	1	2	4	6	12	35

### Table 3: Using a Regulator (0 - 15LPM) with Conventional Cylinders + Using Integrated

Flow setting Continuous Flow	'B' 160L	'CS' 270L	'CH' 470L	Presence PRC TAKE O <sub>2</sub> 590L	'CL' 760L	Presence TAKE O <sub>2</sub> 1000L	'D' 1500L	'E' 4200L	Presence PRE TAKE O <sub>2</sub> 4100L
1	26	4.5	7.8	9	12.6	16	25	70	68
1.5	1.7	3	5.2	6	8.4	11	16.7	46.7	45
2	1.3	2.2	3.9	5	6.3	8	12.5	35	34
3	0.9	1.5	2.6	3.3	4.2	5	8.3	23.3	22
4	0.6	1.1	2	2.5	3.1	4.3	6.2	17.5	17
5	0.5	0.9	1.5	2	2.5	3.2	5	14	13.5
6	0.4	0.7	1.3	1.6	2.1	2.8	4.1	11.6	11.4
9	-	-	-	1	1.4	2	2.8	7.8	7.6
12	-	-	-	0.8	1	1.4	2	5.8	5.7

\* Take 02 is an Oxygen bottle with a built in regulator


# Going out with your Oxygen bottle

When using your Oxygen bottles outside of your home, be aware of your surroundings  
It is your responsibility to check if you need permission to use your Oxygen bottles on public transport and airlines

Using your Oxygen bottle in a vehicle – Only use Oxygen in your car if you absolutely have to.

DO NOT allow smoking in the car even if your Oxygen bottle is turned off

Secure your Oxygen bottle upright in the vehicle so it does not move if there is a sudden stop or accident

Always allow fresh air flow into your car. You can:

- Roll down AT LEAST one window of your car and/or
- Set the vehicle's air conditioning to draw in fresh air

**DO NOT** use your Oxygen when your car is being refuelled

Take only the number of Oxygen bottles you will need for your trip

Take your Oxygen bottles with you when you leave the car. **DO NOT** store Oxygen bottles in your car or the boot

# Fixing problems

## Problem 1: Oxygen is not flowing

Possible cause	Solution
The tubing has come off	Connect the tubing back to the Oxygen bottle
There is a hole in your tubing or tubing is kinked	Un-twist your tubing or replace it as required
The Oxygen bottle is empty	Replace with a full Oxygen bottle
The Oxygen bottle is turned off	Turn your Oxygen bottle on

## Problem 2: Oxygen is leaking and/or you hear a hissing sound

Possible cause	Solution
Regulator is not on properly	<ol style="list-style-type: none"><li>1. Turn off the Oxygen bottle</li><li>2. Take the regulator/OCD off your Oxygen bottle</li><li>3. Put it back on</li><li>4. Check for leaks</li></ol> If the leak and/or hissing does not stop: <ol style="list-style-type: none"><li>5. Close the black 'cap' by turning it to the right</li><li>6. Close the black 'cap' by turning it to the right</li><li>7. Take the regulator/OCD off the Oxygen bottle</li><li>8. Put the faulty regulator/OCD aside</li><li>9. Call Air Liquide Healthcare to get a replacement regulator/OCD</li></ol>
Damaged seal	<ol style="list-style-type: none"><li>1. Turn off flow from the regulator/OCD</li><li>2. Turn off the Oxygen bottle</li><li>3. Remove the regulator/OCD</li><li>4. Check for any visible damage to the seal such as nicks, cuts, discoloration or warping that has built up on it</li><li>5. If the seal is not damaged, put the regulator/OCD back on</li><li>6. Check for leaks again. If the leak persists, discontinue use and call Air Liquide Healthcare</li></ol>
The regulator is loose	Hand tighten the regulator by turning the knob to the right

### Problem 3: OCD does not pulse

Possible cause	Solution
The OCD needs to be reset	<ol style="list-style-type: none"> <li>1. Take the batteries out of the OCD</li> <li>2. Wait 10 seconds</li> <li>3. Put the batteries back into the</li> </ol>
It is not turned on	Press the power button to turn it on
The batteries are flat or weak	Replace the batteries
The batteries are not inserted correctly	<p>Check that the batteries have been inserted correctly</p> <p>The (+) and - on the batteries should line up with + and - on the OCD battery terminals</p>
The strap is blocking the battery holder	Remove the batteries and re-insert them correctly
Oxygen bottle is turned off	Open your Oxygen bottle by turning the black 'cap' on top of the bottle to the left
Oxygen bottle is empty	Replace with a full Oxygen bottle
Tubing is bent or has kinks	Straighten out the tubing and/or replace it if needed
Incorrect batteries are installed	Always replace your batteries with the same model and type of battery originally provided by Air Liquide Healthcare
You are using an Oxygen mask	OCDs can only be used with nasal cannulas
OCD is set to 'continuous flow' mode	Set the OCD to 'Pulse' mode

Call Air Liquide Healthcare if your OCD is not operating properly after troubleshooting . **DO NOT** try and fix it yourself  
**DO NOT** use tape sealer or other material on your OCD to stop leaks

**For further information of the safe use, operation and maintenance of your regulator or OCD, please refer to the supplier's user manual supplied with the device. Supplier user manuals for all devices can also be found on our website [home-oxygen.com.au](http://home-oxygen.com.au) .**

# Paediatric Equipment

---

For paediatric patients, similar types of equipment are supplied as those referenced in the Stationary Oxygen Concentrator and Oxygen Bottle (cylinder) sections, however they will deliver Oxygen flow rates lower than those required or prescribed for adults.

For low flow rates under 2 lpm, a low flow Oxygen concentrator or medical Oxygen cylinder and regulator and/or flow meter may be used to ensure an accurate dosage. The equipment shown below is available:

E cylinder with Nova regulator and flowmeter

E cylinder with Nova regulator and flowmeter

**Always use the prescribed flow rate as per physician instructions.**

# Table 1: Equipment and available flow rates

Paediatric and low flow Oxygen equipment			
Equipment name	ALH Item code	Available flow rates (litres per minute)	Item photo
Mediselect II Paediatric Regulator + any size Oxygen cylinder	6926B	0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8, 1.0, 1.5, 2.0 LPM.	
Alpinox regulator + any size Oxygen cylinder	5481B	0.1 lpm, 0.125 lpm, 0.25 lpm 0.4 lpm, 0.5 lpm, 0.6 lpm, 0.7 lpm 0.8 lpm, 1.0 lpm, 1.5 lpm, 2.0lpm	
DeVilbiss Paediatric 0.125 to 2lpm Litre Oxygen Concentrator	6891B	0.125 lpm, 0.25 lpm, 0.5 lpm 1.0 lpm, 1.5 lpm, 2.0 lpm	

**As with all medical Oxygen, it is very important that the flow rate of Oxygen is administered as per the prescription.**

Generally nasal cannulas are prescribed to be used for paediatric patients, with various sizes available depending on their age. Ensure the following special considerations are followed for paediatric patients:

- Position the nasal prongs along the patient’s cheek and secure the nasal prongs on the patient’s face with adhesive tape.
- Position the tubing over the ears and secure behind the patient’s head. Ensure straps and tubing are away from the patient’s neck to prevent risk of airway obstruction.
- Change the adhesive tape weekly or more frequently as required.

- Check nasal prong and tubing for kinks or twists at any point in the tubing and clear or change prongs if necessary.

**Oxygen masks, bi-flow masks or oxymisers can only be supplied to patients when prescribed by a physician as they can change the dosage of Oxygen.**

Humidification may also be prescribed for pediatric patients to prevent the drying out of the nasal passages. Please see page 28 for information regarding the use of humidifiers.

If at any time you are unsure or have any questions please contact Air Liquide Healthcare on 1300 36 02 02 or call your prescribing physician.

# Oxygen Accessories

## Nasal Cannula

Cannula go in your nostrils and are held in position behind your ears with the tubing secured under your chin. They are:

- Clear and latex free
- Used with all forms of Oxygen therapy
- Suitable for flow rates up to 4 L/min
- Available in a number of different sizes including adult, child, neonate and infant

## Oxygen Masks

These masks fit over your nose and mouth. They are:

- Held into position with an elastic placed around your head
- Clear and latex free
- Are used for flow rates greater than 4L/min
- Available in a number of different sizes including adult, child and infant

## Oxymisers

Oxymisers are Oxygen conserving cannula that use a reservoir system. They need as much as 75% less Oxygen flow to deliver the same amount of Oxygen. They make your cylinders last longer or provide you with a higher percentage of Oxygen at a lower flow rate

- They cannot be used with a humidifier
- We can not give you these without a doctor's prescription

Flow rate in L/min on traditional cannula	Equivalent flow rate in L/min using an oxymiser
2	4
3	5.5
4	6.5
6	8.5


# Oxygen Accessories

Oxygen tubing is available in a variety of different lengths.

1.2M 2.1M 4.6M 9.1M 10.7M

The different lengths are suitable for different types of Oxygen delivery. We can help you choose the right lengths for your Oxygen system.


## Oxygen connectors

Connect your Oxygen tubing to your cannula or mask. They are available as straight or in a 360 degree swivel option.


## Humidifiers

Humidifiers make using Oxygen more comfortable if you have a dry nose or a lot of nosebleeds. They can only be used on stationary Oxygen concentrators and Oxygen bottles.

## Water Traps

Water traps are used to collect condensation inside your Oxygen tubing caused by changes in room temperature.

## Comfort & Protection for Sensitive Skin

EZ Wraps™<sub>1</sub> or Comfy Ears<sub>2</sub> provide a crucial protective barrier for patients experiencing discomfort from medical tubing, devices, or prolonged contact with medical aids.

## NozOil

This special spray is used if the inside of your nose feels dry after you start using Oxygen. Unlike most oils or creams, it is compatible for use with your Oxygen therapy.

"Please note if your Oxygen equipment is being funded by programs such as State Wide Equipment Program (SWEP), Medical Aids and Subsidy Program (MASS), SilverChain, Department of Veterans Affairs (DVA) or Home Care package etc your funding program will include a set quantity of consumables, however for any quantities above the limit these will be need to be purchased privately. Please visit our online store to purchase additional consumables: [store.airliquidehealthcare.com.au](http://store.airliquidehealthcare.com.au)"

# General Oxygen Equipment Cleaning

## Masks, Cannula and Tubing:

### Daily Care:

- **Disconnect:** Turn off the Oxygen and detach the cannula.
- Wipe the nasal prongs with a clean fragrance free wipe daily to remove facial oils.

### Avoid:

Do not use harsh chemicals or alcohol-based cleaners, as they can degrade the tubing.

### Every fortnight:

- **Disconnect:** Turn off the Oxygen and detach the cannula.
- **Replacement:** Replace the cannula every 2 weeks and immediately after respiratory illnesses.
- Change it more often if it become hard, cracked or if your nose gets sore.

### Every 3 months:

- Change your tubing and Oxygen connectors.
- Change your water traps .

## Humidifier:

### Daily Care:

- Turn off the device and disconnect Oxygen source and the power.
- **Empty and Rinse: Remove the humidifier bottle, empty the water, and rinse it with hot water every day.**
- **Wash:** Wash the bottle with mild detergent and rinse with warm water.
- **Water Type:** Always use distilled to avoid mineral buildup, and fill to the recommended level.

### Every 2 months change the humidifier.

### Every week:

1. Turn off the device and disconnect Oxygen source and the power.
2. Separate the parts of your humidifier & wash in warm soapy water.
3. Rinse thoroughly.
4. Soak all parts in a half water, half vinegar mixture for 30 minutes .
5. Rinse thoroughly.
6. **Air Dry in shade:** Do not dry with a towel to prevent contaminating the clean bottle.
7. Put the parts back.
8. **Inspect Regularly:** Check for cracks, discolouration, or wear on the lid and bottle, replacing immediately if found.

## Fire Safe™ Valves:

### Every week:

- **No Cleaning Required:** The valve should not be washed or cleaned internally.
- Should not be submerged or washed like a reusable cannula.
- Take care not to let any water inside the valve .
- **Single Use:** Once a Firesafe valve has stopped a fire (actuated), it cannot be reset and must be discarded.
- **Replacement:** Replace the valve immediately if it becomes contaminated, cracked, or if it is at its expiry date.
- **Alcohol Gels:** If using alcohol hand gel, ensure it is fully evaporated and hands are completely dry before touching the equipment.

# Therapy Installation and Subsequent Visits to a Patient's Home

---

## Using Your Oxygen Safely at Home

### Oxygen Tubing

For your safety and to ensure you receive the correct amount of Oxygen, please keep your tubing length within the following limits:

- Maximum Length: The tubing between your Oxygen machine and your nasal cannula (nose piece) must not exceed 15 meters (49 feet).
- Note: Using longer tubing may drop the pressure and prevent you from getting your prescribed dose.

---

### Where to Place Your Oxygen Equipment

Oxygen is a fire hazard. To keep your home safe, do not install or store Oxygen equipment in the following areas:

- Walkways and hallways: To prevent tripping or blocking exits.
- Kitchens or garages: Stay away from heat sources or where grease and oil are stored.
- Workshops: Avoid areas where sparks may occur (e.g, from power tools).
- Bathrooms: High humidity can damage the equipment.

---

### Safety Records and Paperwork - HHC Risk Assessment

Your Home Healthcare (HHC) Risk Assessment is a vital record of your home's safety setup.

- Record Keeping: We keep a summary of your risk assessments on file for at least five years, or as required by local laws.
- Access: You can request the most recent version of your assessment at any time.

# Therapy Installation and Subsequent Visits to a Patient's Home

---

Our team is dedicated to your ongoing care. We will visit your home for:

1. Initial Installation: Setting up your equipment and teaching you how to use it.
2. Regular Maintenance: Scheduled visits to check your equipment and ensure your therapy is working effectively.
3. Safety Checks: Periodic reviews of your home environment to keep you safe.

## The Customer Service Centre can receive a request to install:

- Equipment for therapy at a new patients home
- Equipment for a new therapy at an existing patients home

## A Home Healthcare (HHC) patient profile is established and includes:

- The physician's name
- The prescription
- The medical device type and setting

## During the first visit to a patient's home, the Home Care Technician (HCT) will:

- Performs, for Oxygen therapy, a technical risk assessment using the checklist in Annex 4.
- Educates the patient, any relatives and/or the caregiver(s) about the therapy.
- Evaluates the effectiveness of the patient's relatives' and/or caregivers knowledge and fill out a form using the template in Annex 5.

## Based on the results of the assessment, the HCT will:

- Either proceeds to the installation
- Decides if corrective action is needed and implements it
- Or refers to the Field Service Manager (FSM) who will determine if corrective action can be implemented or if the installation has to be suspended.

# Therapy Installation and Subsequent Visits to a Patient's Home

---

**The HCT sends the three completed checklists to the FSM who reviews them and establishes a HHC Home Risk assessment summary, including at least:**

- Completed checklists (environmental, technical (if applicable) and patient)
- List of identified issues (if any)
- List of corrective action decided and implemented (if any)
- Final decision to install the equipment or not

**If the HCT detects a change compared with the current HHC risk assessment summary, he/she:**


- Reviews and updates the previously filled checklists (environmental, technical (if applicable) and patient)
- Decides if corrective actions are needed and implements them or refer to the FSM who will determine if corrective action needs to be implemented or if the installation has to be suspended
- Sends the updated checklist to the FSM

**Before a visit to a patients home, the HCT reviews the last HHC risk assessment summary. During a visit, the HCT:**



- Determines if there have been any changes since the previous visit by reviewing the risks related to the therapy and updates the assessment if necessary:
  - The environment risk assessment using the checklist in Annex 3;
  - The technical risk assessment for Oxygen therapy using the checklist in Annex 4;
- Evaluates the patient's, any relatives' and/or caregiver's knowledge and/or understanding of the therapy using the checklist in Annex 5;
- If needed, refreshes the patient's, relatives', and/or caregiver's knowledge.

**The HCT reviews the updated checklists and if needed updates the HHC risk assessment summary that includes at least:**



- The completed checklists (environmental, technical (if applicable) and patient);
- The list of identified issues (if any);
- The list of corrective action decided and implemented (if any);
- The final decision to install the equipment or not.

# O2 Oxygen Storage Safety

	<p>The cylinders are stored at a minimum distance of 3 m (10 ft) away from open flames.</p>
	<p>The cylinders are stored at a minimum distance of 1.5 m (5ft) away from electrical equipment.</p>
	<p>The cylinders are stored at a minimum distance of 1.5 m (5 ft) away from combustible material.</p>
	<p>Oxygen cylinders stored outside the patient's home are located in a well-ventilated, safe and covered area and not in a closed space or cupboards</p>
	<p>Two storage areas are identified, one for empty cylinders) and one for full cylinders.</p>
	<p>Oxygen cylinders stored inside the patient's home are installed near an exit and not blocking it.</p>
	<p>Oxygen cylinders stored inside the patient's home are secured to prevent them from falling.</p>
	<p>Portable cylinders used for patient mobility are stored on their side or in a specifically designed rack.</p>

# O2 Oxygen Storage Safety

	<p>The device is installed at a minimum distance of 3 m (10 ft) away from open flames.</p>
	<p>The device is installed at a minimum distance of 1.5 m (5 ft) away from electrical equipment.</p>
	<p>The device is installed at a minimum distance of 1.5 m (5 ft) meters away from combustible material.</p>
	<p>The device is installed at a minimum distance of 0.15 m (0.5 ft) meters away from curtains, upholstered seats or any other material that can prevent adequate air circulation around the device.</p>
	<p>The device is installed in such a way that the concentrator's air intake or exhaust is not obstructed.</p>
	<p>The device is installed close enough to the patient for them to hear the alarms</p>
	<p>The device is installed close to a wall electrical outlet.</p>
	<p>The device is installed away from water outlets and away from areas where water might be spilled or could condense.</p>
	<p>The device is installed in a room with a low or normal humidity level.</p>
	<p>A self-filling system is installed with the same rules as for the installation of a standard Oxygen concentrator in a place where it can be operated safely</p>

# Oxygen Reference Guide

## General Safety Precautions

### 1. Always read instructions

Read the User **Instruction Manual** carefully before operating your oxygen equipment. Pay special attention to information where hazard symbol are shown.

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### 2. No smoking. No naked flames.

Materials burn much more vigorously in oxygen than in air. Never smoke (or let someone else smoke near you) whilst using your oxygen equipment. Do not use your oxygen equipment within 3 metres of open fires or naked flames. Post at least one 'No Smoking - Oxygen in Use' sign in a prominent place.

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### 3. Use in a well ventilated area

Only use and store your oxygen equipment in a well ventilated area. Keep internal doors open whilst your medical oxygen is in use.

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### 4. Avoid oxygen enrichment

Never place your oxygen equipment near curtains or cover it with clothing, as this will restrict air circulation. Never use or carry your oxygen equipment under any clothing or bedding. Fire hazards are increased if materials become oxygen enriched with no ventilation.

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# General Safety Precautions

## 5. Secure and use oxygen equipment in an upright position

Follow the instructions provided by Air Liquide Healthcare for the safe storage and use of oxygen equipment. Oxygen cylinders must be secured upright to prevent them from falling over. Oxygen concentrators must always be stored and transported upright.

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## 6. Never use oil or grease

Do not use oil or grease with your oxygen equipment. Only use skin moisturizer products that are free of petroleum based ingredients.

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
## 7. Oxygen tubing

Attach the oxygen tubing to the outlet of your oxygen delivery system or FireSafe™ Valve. For equipment delivering continuous flow oxygen, the length of tubing should not exceed 15 metres. Most oxygen conserving devices and portable oxygen concentrators only use a maximum tubing length of 2 metres. Always verify the maximum tubing length that can be used in the applicable manufacturers' 'User Instruction Manual'.

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## 8. Adjust flowrate as prescribed

Set the flow control valve to the flow rate **prescribed by your Doctor**. For accuracy, read the flow rate at eye level.

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## 9. Phone for assistance if required

If for any reason your oxygen equipment fails, consult the troubleshooting guide. If further assistance is required call Air Liquide Healthcare immediately. Never try and repair any fault unless specifically instructed by Air Liquide Healthcare.

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# Oxygen Concentrator Safety Precautions

## 10. Turn on to obtain oxygen

Plug your oxygen concentrator into the electrical wall outlet supply and switch on. Never connect your concentrator to double adaptors or power boards.

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
## 11. Using FireSafe™ Valves

FireSafe Valves are used to stop the flow of oxygen in the event of fire. Use two FireSafe valves: one at the outlet of the stationary concentrator and the second near you to connect your oxygen cannula /mask to the extension tubing. Some models of concentrators have built in FireSafe Valves and you will only need to use the one valve nearest you. Your technician will advise you if you need one or two FireSafe valves.

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## 12. Keep filter clean

**Regularly** clean the air inlet filter as instructed in the specific manufacturer 'User Instruction Manual'.

Ensure that the filter is **dry** before reusing. **Switch off** and isolate your oxygen concentrator from the mains supply when cleaning or replacing the filter.

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## 13. Only clean with a cloth

Follow the 'User Instruction Manual' for regular cleaning of your oxygen concentrator and associated accessories. Only wipe the exterior with a clean damp cloth (no abrasives) and allow the equipment to dry completely before using.

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15. Do not open the concentrator

**Never** open or remove the concentrator cover at any time. **Only** qualified Air Liquide Healthcare personnel are authorized to carry out repairs on your oxygen concentrator. Doing otherwise may cause electrocution.

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## Oxygen Cylinder Safety Precautions

16. Prepare the Cylinder

Place the cylinder in a trolley or carry bag. Remove the tamper-evident seal and the Valve plug from the cylinder valve. **Slowly** open the cylinder valve counter clockwise to blow out any dust that may have accumulated, ensuring the valve outlet opening is pointing away from you and from other persons when doing so.

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# Oxygen Cylinders Quick Reference Guide

## 17. Do not use excessive force

Place the regulator or oxygen conserving device over the cylinder valve. Align the pins on the regulator or oxygen conserving device with the holes on the cylinder valve and hand tighten in position. **Never** use excessive force or tools as over tightening may damage the equipment.

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
## 18. Do not use the cylinder with a regulator/OCD attached if a leak persists

Open cylinder valve **slowly** by turning it counter clockwise as far as it will go, then back a quarter turn. If leaks occur between the regulator/oxygen conserving device and cylinder, hand tighten until the leak stops. If the leak persists, discontinue the use of the equipment and contact Air Liquide Healthcare immediately.

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
## 19. Check cylinder contents

Regularly check how much oxygen is available by checking the gauge on the regulator or oxygen conserving device. When the contents gauge pointer is in the red zone, exchange your cylinder in use for a full one as soon as possible.

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# POWER FAILURE

## EMERGENCY PLANNING

Air Liquide Healthcare's primary concern is for your safety. If you have lung disease, you may be at risk during an emergency or disaster if there is an electrical power failure and your Oxygen concentrator cannot be used.

You are encouraged to develop an emergency plan in the event of power failure. Everyone in your family, including your physician, carers and neighbours should be included in this plan. It is recommended that you

review your Oxygen emergency plan together every 6 months or more frequently if there is a change in your Oxygen needs.

You should document your emergency Oxygen plan and keep it readily available so that you can refer to it when needed. Your relatives / carers / neighbours / physician / friends may assist you in developing this emergency plan.

### PREPARING FOR AN EMERGENCY

Use this checklist to help you prepare your Oxygen emergency plan of action in the event of a power failure

- |   |  |
|---|--|
| <input checked="" type="checkbox"/> My electricity provider has been advised that I am dependent on electricity to supply my Oxygen.    | <input checked="" type="checkbox"/> I have the names and contact numbers of my relatives/ carers/neighbours/physician/ friends up to date and documented.  |
| <input checked="" type="checkbox"/> I check on a weekly basis that my emergency Oxygen back-up cylinder is full.                        | <input checked="" type="checkbox"/> I have the names and telephone numbers of support services I may require written on my plan. These may include: police, state emergency service, Ambulance, community Health centre, doctor(s), family, Air Liquide Healthcare and anyone else you deem important that can assist. |
| <input checked="" type="checkbox"/> I have torches with spare batteries available for emergency use.                                    | <input checked="" type="checkbox"/> I will date my plan and review it every 6  |
| <input checked="" type="checkbox"/> If a telephone is not available, alternate assistance has been arranged.                            | <input checked="" type="checkbox"/> This plan is easily accessible in case of emergency.   |
| <input checked="" type="checkbox"/> I have a radio with spare batteries available for emergency use to hear electricity supply updates. |  |

### IMPORTANT INFORMATION TO NOTE

Please advise Air Liquide Healthcare of the following points:

- If your home address changes (even temporary change) so our Home Healthcare Technicians are aware, they may not be delivering to your main address.
- If your Oxygen therapy has been interrupted or if you switch to treatment.
- If you are admitted to the hospital.

# MY EMERGENCY PLAN

## My power failure Oxygen emergency plan:

**Step 1:** use the torch if lighting is required - Never use matches or candles.

**Step 2:** shut the power off at the outlet of the stationary Oxygen concentrator.

**Step 3:** use the back-up or portable Oxygen cylinder.

**Step 4:** select the flow rate on the regulator connected to the Oxygen cylinder to \_\_\_\_\_ Lpm.

**Step 5:** contact the electricity supplier to see how long the power failure will last.

**Step 6:** contact the appropriate emergency contact/s.

**Step 7:** Limit my physical activity and stay calm and relaxed.

**Step 8:** Listen to the radio for updates.

## My emergency contacts

CONTACT	COMPANY/NAME	PHONE NO.
Police		000
State emergency services		
Ambulance		000
Community Health Centre		
Electrical supplier		
Primary physician		
Primary carer		
Neighbour		
Friend		
Home Oxygen supplier	Air Liquide Healthcare	1300 36 02 02

If a telephone is unavailable my plan is: \_\_\_\_\_

Other important information: \_\_\_\_\_

Changing Care with  
Australia since 1995

**70+**  
Clinic &  
Partner sites  
Nationally

**60,000**  
Patients  
treated in  
their Home

**+ 200**  
Employees

**+100,000**  
OSA Patients  
Supported

## CONVENIENTLY LOCATED ACROSS AUSTRALIA

For more information and product manuals  
visit our website at

<https://au.healthcare.airliquide.com/home-oxygen-products>

### HEAD OFFICE

Suite 4/270 Lahrs Road  
Ormeau, QLD 4208