

# Relative Analgesia



## Nitrous oxide sedation equipment for dental practices

Specifically designed nitrous oxide sedation equipment and administration systems

Air Liquide Healthcare offers a range of equipment and administration systems specifically designed for dental surgery set ups.

- ✓ Range of options catering for a variety of needs
- ✓ Solutions compliant with industry recommendations

# Nitrous oxide / oxygen mixers

## Fraser MDM®

- Mechanical mixer
- Auto-compensation system
- Oxygen fail-safe and flush

Part # 583S



## Digital MDM®

- Digital mixer
- Audible and visual alarms
- Auto-compensation system
- Oxygen fail-safe and flush

Part # 3360S



## Mounting options

### 2-cylinder mobile trolley

- Accommodates 2 D-size cylinders
- Provides complete mobility
- Ideal for surgeries without piped gas



H: 94 cm  
W: 50 cm  
Depth: 35 cm  
Part # 592S

### Telescopic pedestal

- Adjustable height
- Mobile solution for surgeries equipped with piped gas



H(min): 70cm  
H(max): 98cm  
Radius: 65cm  
Part # 3335S

### Wall and cupboard mounting

- Special "Slide-out" design for cupboards
- Suitable for surgeries equipped with piped gas



15" telescopic arm  
Part # 605S



8" arm  
Part # 4809S

Cupboard Slide Assembly  
Part # 4810S

## Patients circuits and nasal hoods

### Matrx Scavenger

- Complete scavenging patient circuit
- Fully autoclavable
- Comes with a large nasal hood



Part # 3570S

### Reusable nasal hoods

- Suitable for Matrx Scavenger
- Fully autoclavable
- 3 sizes available



Small: Part # 848S  
Medium: Part # 847S  
Large: Part # 846S

### Disposable scented nasal hoods Dynomite™

- Suitable for Matrx Scavenger
- 5 fresh fragrances
- 3 sizes available



Small: Part # 1546S  
Medium: Part # 1547S  
Large: Part # 1548S  
Assorted: Part # 1549S

## Other Equipment

### Breathing Bag

- 3 litre capacity
- Latex-free



Part # 4882S

### Regulators



Medical oxygen regulator  
Part # 4542S



Medical nitrous oxide regulator  
Part # 944S

### Hand wheel and hose assemblies



Medical oxygen hose  
Part # 4737S



Medical nitrous oxide hose  
Part # 4736S

### Mounting pin

For use with 2 cylinder mobile trolleys and wall mounting options



3" pin  
Part # 1522S