

# Which Portable Oxygen Concentrator (POC) is the best for *you?*



| Model                                                                                                    | Type of POC   | Modes of O2                   | Flow Setting                                                    | Maximum Oxygen Output^                                | Weight* | POC run time (at setting of 2*)                                                                                                                                  | Level (dBA)                      |
|----------------------------------------------------------------------------------------------------------|---------------|-------------------------------|-----------------------------------------------------------------|-------------------------------------------------------|---------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------|
| Simply Go Mini          | Portable      | Pulse flow                    | Settings 1-5 pulse flow                                         | 1 LPM                                                 | 3.7 kg  | Comes standard with 1x extended battery<br>4.5 hours (standard battery)<br>9 hours (extended battery)<br>Best for: Daily chores and family visits                | 43*                              |
| Rove 6 16-Cell          | Portable      | Pulse flow                    | Settings 1-6 pulse flow                                         | 1.2 LPM                                               | 3.3 kg  | Comes standard with 1x extended battery<br>5 hours (standard battery)<br>10 hours (extended battery)<br>Best for: Active lifestyle and travel                    | 39*                              |
| Inogen® One G5 16-cell  | Portable      | Pulse flow                    | Settings 1-6 Pulse flow                                         | 1.2 LPM                                               | 3.3 kg  | Comes standard with 1x extended battery<br>5 hours (standard battery)<br>10 hours (extended battery)<br>Best for: Active lifestyle and travel                    | 38**                             |
| Simply Go              | Transportable | Pulse flow<br>Continuous flow | Settings 1-6 pulse flow<br>Settings 0.5 - 2 lpm continuous flow | 1.4 LPM (on pulse dose)<br>2 LPM (on continuous flow) | 8.2 kg  | Comes standard with 1x rechargeable battery<br>Pulse Flow: 3.5 hrs. Continuous Flow: 0.7 hours<br>Best for: Single unit solution that for all-around daytime use | 43* (pulsed)                     |
| Eclipse 5             | Transportable | Pulse flow<br>Continuous flow | Settings 1-9 pulse flow<br>Settings 0.5 - 3 lpm continuous flow | 2.3 LPM (on pulse dose)<br>3 LPM (on continuous flow) | 11.5 kg | Comes standard with 2x rechargeable batteries<br>Pulse flow: 5 hours<br>Continuous flow: 1 hour 48 min<br>Best for: Those needing higher flow rates of oxygen    | 40* (pulsed)<br>48* (continuous) |
| Zen-O                 | Transportable | Pulse flow<br>Continuous flow | Settings 1-6 pulse flow<br>Settings 0.5 - 2 lpm continuous flow | 2 LPM (on pulse dose)<br>2 LPM (on continuous flow)   | 6 kg    | Comes standard with 2x 12-Cell batteries<br>Pulse flow: 4 hours (single battery)<br>Continuous flow: Up to 45 min (single battery)<br>Best for: Active patients  | 42*                              |

\*POC run time is based on in-house testing performed by Air Liquide Healthcare at a setting of 2, and is based on 18-20 breaths per minute. Duration times of optional external batteries are based on manufacturer specifications. \*Weight inc bag/cart/accessories. \*Sound level based on in-house testing performed by Air Liquide Healthcare. \*\*Sound level based on manufacturer specifications; ^calculated at 20 breaths per min at setting of 2.

Always read the instructions for use before operating a medical device.

Please do not hesitate to contact us for any assistance. We will help you find the POC that meets your oxygen needs and complements your way of life.

 [au.healthcare.airliquide.com/vital-o2](http://au.healthcare.airliquide.com/vital-o2)  [store.airliquidehealthcare.com.au](http://store.airliquidehealthcare.com.au)  1300 360 202  [alhenquiries@airliquide.com](mailto:alhenquiries@airliquide.com)

# Making Portable Oxygen *easy*

Whether you're purchasing or renting a portable oxygen concentrator (POC) from Air Liquide Healthcare, we ensure peace of mind with every step.



## **Purchase a portable** *Oxygen concentrator*

Buying a POC may seem like a big decision, but we're here to provide all the necessary information for a smooth process. Enjoy peace of mind with the following:

- ✓ Complete training on safe use<sup>^</sup>
- ✓ Equipment to meet all Medical Oxygen Therapy needs
- ✓ Compliant with Medical Standards
- ✓ Wide selection of accessories
- ✓ Dedicated customer service and aftercare
- ✓ Qualified Technicians for all servicing needs
- ✓ Trusted by Physicians and Health Care Providers
- ✓ Airline Approved Equipment
- ✓ Finance options available

## **Rent a portable** *Oxygen concentrator*

Short-term or long term rental oxygen needs? We've got you covered!

If you require a Portable Oxygen Concentrator for rental, we offer a range of rental options, to get started on one of our rental plans:

- ✓ Quick and easy rental options
- ✓ Long-term or short term rentals available\*
- ✓ Full support from our team

### **Get started today;**

- 1 Speak with your medical oxygen provider to determine the right solution and have a prescription complete
- 2 Obtain our medical prescription online at <http://au.healthcare.airliquide.com/oxygen>

Once completed please fax your form to

**1800 254 329** or email to:

**[alorders@airliquide.com](mailto:alorders@airliquide.com)**

Our team will contact you to discuss your options.

You may rent a specific model to try before you purchase upon availability, min. days apply. If you decide to proceed with a purchase post trial of the same model, you will be eligible for a discount off the purchase price. <sup>^</sup>Home Care Technician In-Home Set Up & Patient Education support at additional cost. \*Short Term Rental POCs available from minimum of 2 wks. Long Term Rental POCs available from 6 mths.

**Always read the instructions for use before operating a medical device.**